



Animal Health Laboratory Services
Hancock Building – 65 River Rd, Truro, Nova Scotia, B2N 2P3

MILK ANTIBIOTIC RESIDUE TESTING

Client/Producer Information	
Name	
Address	
Phone/Fax	
Email	
Report to	<input type="checkbox"/> Client <input type="checkbox"/> Other: _____
Species Information: <input type="checkbox"/> Cow <input type="checkbox"/> Goat/Sheep	
Collection Time and Date:	
No.	Sample#/ID
#1	
#2	
#3	
#4	
#5	
#6	
Lab Use Only:	
VADDS Entry by:	
Received in Lab by:	

LAB USE ONLY		
AHL Laboratory ID	Case Coordinator	
Charm Trio (Cow) <input type="checkbox"/>	SL(BL) Goat/Sheep <input type="checkbox"/>	
Test Strip lot#:		
Analyst:	Date tested:	
Calibration and Set up Information		
Low Calibrator reading (-42 to -28)		
High Calibrator reading (152 to 226)		
Incubator Temperature (56 ± 1 °C)		
Thermometer SN:	°C	
Expiry:		
Pipette Verification (0.300g ± 0.015g)		
Pipette used <input type="checkbox"/> 023225 <input type="checkbox"/> 077634	g	
Balance ID: _____		
Controls		Results
Batch ID: _____	Negative < -600	
	Positive > +400	
Sample#/ID	Results	Drug Identified
Client notification		
Telephone notification		
Date/initials		
Verification and Reporting		
Date/Initials		

Specimens are analyzed as provided. The laboratory takes no responsibility for the accuracy of the information.